MEDISANA®

## Elektrischer Nagelpolierer NP 860

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# MEDISANA®

## Electric nail polisher NP 860

- ① Battery compartment
- ② Slide switch
- ③ Attachment

### Explanation of symbols

This instruction manual belongs to this device. It contains important information about starting up and operation. Read the instruction manual thoroughly. Nonobservance of these instructions can result in serious injury or damage to the device.

#### WARNING

These warning notes must be observed to prevent any injury to the user.

#### ATTENTION

These notes must be observed to prevent any damage to the device.

#### NOTE

These notes give you useful additional information on the installation or operation.

**LOT** LOT number

Manufacturer

## GB Safety Instructions

**Read the instruction manual carefully before using this device, especially the safety instructions and keep the instruction manual for future use. Should you give this device to another person, it is vital that you also pass on these instructions for use.**

- Only use this device for the purpose for which it was developed and in the manner specified in the instructions for use. When using it for other purposes, the warranty expires.
- This device is not designed for commercial use.
- Do not use any auxiliary parts that have not been recommended by the manufacturer.
- This appliance can be used by children aged from 8 years and above and persons with reduced physical, sensory or mental capabilities or lack of experience and knowledge if they have been given supervision or instruction concerning use of the appliance in a safe way and understand the hazards involved.
- Children must be supervised to ensure that they do not play with the device.
- Do not attempt to repair the unit yourself in the event of any defects. Not only does this render the warranty null and void but you can also run the risk of serious hazards (fire, electric shock, injury). Only have repairs carried out by authorised service points.
- Only use the unit indoors.
- Do not use any electrical equipment while in the bath or shower. Do not store or keep any electrical equipment in places where it could fall into the bathtub or washbasin.
- Keep the device away from sources of heat.
- Never use the device under a blanket, pillow or cushion, etc.
- Uninterrupted treatment should not last longer than 15 minutes. Then switch off the unit for at least half an hour and allow it to cool down.
- Consult your doctor first before using this unit if you have diabetes or other related illnesses.
- The treatment should be pleasant. Stop at once and consult your doctor if you feel pain or the treatment seems unpleasant.
- For hygiene reasons, the attachments should only be used by one person.
- Cleaning and user maintenance must not be performed by children unless supervised.

### BATTERY SAFETY INFORMATION

- Do not disassemble batteries!
- Clean the battery and device contacts if necessary before inserting the batteries!
- Remove discharged batteries from the device immediately!
- Increased risk of leakage, avoid contact with skin, eyes and mucous membranes! If battery acid comes in contact with any of these parts, rinse the affected area with copious amounts of fresh water and seek medical attention immediately!

- If a battery has been swallowed seek medical attention immediately!
- Insert the battery correctly, observing the polarity!
- Keep the battery compartment well sealed!
- Remove the battery from the device if it is not going to be used for an extended period!
- Keep batteries out of children's reach!
- Do not attempt to recharge batteries! **There is a danger of explosion!**
- Do not short circuit! **There is a danger of explosion!**
- Do not throw into a fire! **There is a danger of explosion!**
- Keep unused batteries in their packaging away from metal objects in order to prevent short circuiting!
- Do not throw used batteries into the household refuse; put them in a hazardous waste container or take them to a battery collection point, at the shop where they were purchased!

### Items supplied

Please check first of all that the device is complete and is not damaged in any way. In case of doubt, do not use it and contact your dealer or your service centre.

The following parts are included:

- 1 **MEDISANA** Electric nail polisher **NP 860**
- 4 attachment parts (each 2 smoothing and 2 polishing rolls)
- 2 batteries AA LR6, 2x1,5V
- 1 instruction manual

The packaging can be reused or recycled. Please dispose properly of any packaging material no longer required.

If you notice any transport damage during unpacking, please contact your dealer without delay.

#### WARNING

**Please ensure that polythene packing is kept away from the reach of children! Risk of suffocation!**

### GETTING STARTED - Removing / inserting batteries:

Open the battery compartment ① and insert the two supplied batteries of size AA (LR6) 1.5V. Pay attention for correct polarity according to the marking. Do not use rechargeable batteries. Close the battery compartment ①. If the battery power gets low after some time of usage, you need to replace the batteries.

### USE

This device is intended only for the treatment of toenails (pedicure) or fingernails (manicure). Put the desired attachment (smoothing roll = coarse pink roll or polishing roll = smooth, white roll) onto the device. Slide the slide switch ② to first position to run the device with normal speed or to second position to run it with high speed.

Slide the slide switch ② back to its initial position to switch the device off.

**Smooth nails:** Gently smooth the nail surface by slightly pressing the smoothing roll onto it for approximately 15 seconds.

**Polish nails:** Polish the nail surface using light pressure with the polishing roll.

### Cleaning and maintenance

- Clean the unit and the attachments with a soft cloth which has been lightly moistened with a mild soap water or with plain water. Do not let water get into the unit. Only use the unit again when completely dry.
- Never use strong cleaning agents or stiff brushes.

### Disposal

This product must not be disposed of together with domestic waste. All users are obliged to hand in all electrical or electronic devices, regardless of whether or not they contain toxic substances, at a municipal or commercial collection point so that they can be disposed of in an environmentally acceptable manner. Please remove the batteries before disposing of the device / unit. Consult your municipal authority or your supplier for information about disposal.

### Technical specifications

Name and model:	<b>MEDISANA</b> Electric nail polisher <b>NP 860</b>
Power supply:	2x batteries AA LR6, 2x1,5V
Rounds per minute:	approx. 2,000 (level 1); approx. 2,400 (level 2)
Dimensions:	approx. 15.5 x 3.8 x 3.5 cm
Weight:	approx. 150 g incl. batteries
Article No.:	85160
EAN No.:	4015588 85160 5
Spare parts:	4 pcs. pack ± 2 x smoothing roll & 2 x polishing roll; Art no.: 85164; EAN: 4015588 85164 3

**In accordance with our policy of continual product improvement, we reserve the right to make technical and optical changes without notice.**

The current version of this instruction manual can be found under [www.medisana.com](http://www.medisana.com)

### Warranty and repair terms

Please contact your dealer or the service centre in case of a claim under the warranty. If you have to return the unit, please enclose a copy of your receipt and state what the defect is. The following warranty terms apply:

1. The warranty period for **MEDISANA** products is three years from date of purchase. In case of a warranty claim, the date of purchase has to be proven by means of the sales receipt or invoice.
2. Defects in material or workmanship will be removed free of charge within the warranty period.
3. Repairs under warranty do not extend the warranty period either for the unit or for the replacement parts.
4. The following is excluded under the warranty:
  - a. All damage which has arisen due to improper treatment, e.g. non-observance of the user instructions.
  - b. All damage which is due to repairs or tampering by the customer or unauthorised third parties.
  - c. Damage which has arisen during transport from the manufacturer to the consumer or during transport to the service centre.
  - d. Accessories which are subject to normal wear and tear.
5. Liability for direct or indirect consequential losses caused by the unit are excluded even if the damage to the unit is accepted as a warranty claim.

**MEDISANA AG**  
Jagenbergstraße 19  
41468 NEUSS  
Germany  
E-Mail: [info@medisana.de](mailto:info@medisana.de)  
Internet: [www.medisana.com](http://www.medisana.com)

The service centre address is shown on the attached leaflet.

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## INSTRUCTIONS DE SECURITE RELATIVES AUX PILES

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